

---

**Table S5. Primers and adapters used in the study**

	<b>primer name</b>	<b>5'-primer sequence-3'</b>
Inverse PCR Primers for pFind1 colonies 1st Round	PB36inv5F	CAGCGACACCCACCGCACCGACAAGCACGC
	PB152	CCGCAAGCCACGCAACGCGGAACCC
Inverse PCR Primers for pFind1 colonies 2nd Round	PB149	GCGACCGAGACGCCCCAAACGCAC
	PB153	GGGCCACGAACCAACGACCCCGCAA
Ion Torrent adapters for next generation sequencing	A adapter	CCATCTCATCCCTGCGTGTCTCCGACTCAG
	P1 adapter	CCTCTCTATGGGCAGTCGGTGAT